



**Department of Electronic & Telecommunication
Engineering**

University of Moratuwa

**BM2210 - Biomedical Device Design
Final Report**

Group Members

- 230218G GUNATHUNGA U.A.
- 230487D PERERA W.A.L.S.
- 230012U ABEYWARDHANA T.C.W.
- 230481E PERERA K.W.A.O.V.

November 30, 2025

Contents

1	Introduction	2
1.1	Need Statement	2
1.2	Problem Status	2
2	Concept Generation and Selection	3
2.1	Initial Concept Selection	3
2.2	Final Concept and Proposed Solution	4
3	Implementation Methodology of the Proposed Solution	4
3.1	Operating Mechanism of the IR-Based Drop Sensing	4
3.2	Functional Block Diagram	5
3.3	Component Selection	5
3.4	3D Enclosure Design	6
3.4.1	Initial Sketch	6
3.4.2	Evolution of the 3D Design	7
3.5	IoT-Enabled Web Platform for IV Fluid Monitoring	8
3.5.1	Backend	8
3.5.2	Frontend	8
3.6	User Interface	9
3.7	Final Product and Product Dissection	14
4	Accuracy Testing Results	16
5	Final Budget and Pricing Strategy	17
6	Possible Regulatory Pathway	18
7	Conclusions and Future Improvements	19
7.0.1	Conclusion	19
7.0.2	Future Improvements	19
8	Task Allocation	21
9	Program for ESP32 Microcontroller	21

1. Introduction

1.1. Need Statement

Key Problems Identified

1. Flow rates are manually set; deviations cause overdose/underdose risks.
2. Educated patients could help monitor flow if external displays exist.
3. Some patients self-adjust flows (dangerous for chemotherapy).
4. Regulators are expensive and rare in Sri Lanka.

Population

- Medical professionals— doctors, nurses, anesthetists, emergency staff.
- Patients— cancer, post-surgery, ICU, dehydrated, comatose/unresponsive patients

Outcome

Patients and medical professionals in the ward can adjust the IV fluid drip rate to the original set value by identifying possible deviations as indicated in the display.

Final Need Statement

There is a critical need for a safe, reliable, affordable, cost effective, doctor-friendly and patient-friendly method to externally regulate and display IV fluid drip rate in hospitalized patients.

1.2. Problem Status

The problem status can be explained with a user scenario as follows.

In a busy government hospital ward, a nurse or ward master is responsible for administering IV fluids to multiple patients at the same time. With traditional gravity-based infusion pumps, the ward master must manually set a timer and adjust the drip rate—often in increments of 3 drops per minute. However, patient movement or other disturbances can cause the flow to deviate, and since these pumps lack an external display, the ward master must physically check each patient every two hours to ensure the rate is still correct. This constant vigilance is time-consuming and stressful, and even then, deviations may go unnoticed.

In some cases, patients intentionally increase the drip rate—such as chemotherapy patients wanting treatment to finish early—creating further risk.

With our solution, the infusion pump displays the current drip rate in real time on an external screen. If the rate deviates from the user-set target, the system immediately alerts the ward master with an alarm. Even patients who try to adjust the rate themselves cannot exceed safe limits without triggering an alert. The ward master can monitor all patients' infusion rates simultaneously through a mobile app, reducing the need for constant physical checks. Additionally, if the rate drifts, the patient or nurse can manually adjust it using the control wheel while observing the display, ensuring the correct flow at all times.


This approach improves patient safety, reduces the workload and stress on hospital staff, and brings the benefits of advanced infusion pumps to a system that is affordable and practical for government hospitals.

2. Concept Generation and Selection

2.1. Initial Concept Selection

We selected three initial concepts to move forward with:

1. Smart Gravity-Based Infusion Pump with IR-Based Drop Sensing
2. Smart Gravity-Based Infusion Pump with Capacitive Drop Sensing
3. Smart Gravity-Based Infusion Pump with Image Processing-Based Drop Sensing

Concept	Introduction	Testing Results	Advantages	Disadvantages
IR-Based Drop Sensing	Basic principle of detecting droplets using infrared beam interruption.	Initial tests show consistent detection for clear and transparent liquids with a sensitivity of 1 drop per 200 ms. 	Simple setup, low cost, fast response time.	Affected by ambient light (sensitivity can be adjusted to avoid)
Capacitive Drop Sensing	Uses change in capacitance caused by presence of a droplet.	Initial testing showed that although the sensor responds to capacitance changes, obtaining clear and repeatable readings requires further calibration. The system also cannot reliably detect the exact moment a water droplet makes contact, and because the droplet volume is very small, the resulting capacitance change is minimal. Additional tuning and signal conditioning are therefore needed for accurate droplet-arrival detection.	High sensitivity for higher drop rates, works with transparent fluids.	Requires calibration, performance varies with humidity and temperature.
Image Processing-Based Drop Sensing	Detects droplets using camera images and computer vision algorithms.	Achieved high accuracy and detailed analysis in test runs.	Very accurate, provides size/shape information, flexible for many setups.	High computational cost, expensive hardware.

2.2. Final Concept and Proposed Solution

Final Concept: Smart Gravity-Based Infusion Pump with IR-Based Drop Sensing

Proposed Solution: In our proposed solution, we aim to make infusion therapy more accessible and affordable by enhancing a simple, gravity-based infusion pump with the essential features listed below. Existing solutions are often inadequate for government hospitals, including specialized centers like cancer hospitals, due to the lack of external rate displays and monitoring—even in the most advanced infusion pumps currently used in the hospital system. Additionally, the high cost of current devices makes them unattainable for many healthcare facilities.


1. Mechanism to externally adjust the drip rate with real-time indication of drip rate in drops/min and mL (doctor-friendly and patient-friendly)
2. External display of drip rate in drops/min and mL
3. Alarm system for rate deviations
4. Mobile app/IoT-based monitoring system for real-time tracking across multiple beds in a ward
5. Mobile app/IoT-based monitoring system for real-time adjusting of the drip rate across multiple beds in a ward (for example: the ward master can remotely adjust the drip rate by identifying deviations from the originally set value)

3. Implementation Methodology of the Proposed Solution

3.1. Operating Mechanism of the IR-Based Drop Sensing

In this method, a standard infrared (IR) sensor module is used and rewired, so that the IR emitter and receiver are placed directly opposite each other across the IV drip chamber just before the fluid enters the tube. When no drop is present, the IR beam passes uninterrupted and is detected continuously. Each falling drop interrupts the beam, producing a momentary signal change that is captured by the photodetector. These interruptions are counted by an ESP32 microcontroller, which converts the drop count into clinically relevant infusion rates (drops/min) and displays on a seven-segment display. This design requires only three core components, ESP32 module, IR module, and OLED screen that makes this method simple, low-cost, and highly reliable and ideal for affordable hospital use.

3.2. Functional Block Diagram




3.3. Component Selection

Component	Role in the Product
ESP32	The brain of the system. Processes signals from the IR sensor, calculates drop rate, displays results, and connects with mobile/IoT systems.
Buzzer	Acts as an alarm. Alerts staff immediately when the drip rate deviates from the expected value.
Num Pad	Allows medical staff to input the desired drip rate (drops/min) for each patient.
Display	Shows the real-time drip rate in drops/min and mL, making it doctor-friendly and patient-friendly.
TP4056	Manages battery charging with built-in overcharge protection, ensuring safe and reliable backup power.
DC Jack	Provides AC power connection for continuous operation in the ward.
IR Sensor	Detects each IV fluid drop when it interrupts the light beam, enabling accurate real-time drop counting.

3.4. 3D Enclosure Design

3.4.1 Initial Sketch



3.4.2 Evolution of the 3D Design

		
Overall Dimensions	Larger footprint, occupies more volume	Significantly smaller footprint, optimized volume
Internal Component Arrangement	Relaxed internal spacing, easy routing of wiring and PCBs	Tighter internal layout, requires careful component placement
Portability	Relatively heavy and bulky, not convenient for frequent relocation	Lightweight and highly portable, easier bedside handling
Ease of Assembly & Maintenance	Good component access, straightforward assembly and servicing	Moderate compact layout reduces accessibility and requires planned service access points
Material Usage	Higher material consumption (larger panels and supports)	Reduced material usage leading to lower per-unit material cost
Aesthetic & Ergonomics	Boxy and functional, less bedside friendly aesthetic	Sleeker, modern appearance, better suited to clinical bedside environments
Mounting Options	Requires a larger mounting surface or custom bracket	Easier to mount on IV poles, bed frames, or compact bedside brackets
Stability & Rigidity	Inherently stable due to larger base and mass	Needs optimized structural ribs/walls to retain rigidity in a smaller form factor
Clinical Usability	Buttons and labels can be spaced for clarity and to prevent accidental presses	Controls are compacted layout must avoid accidental activation and remain legible
Battery & Electronics Capacity	Room for a larger battery pack and extra modules for expansion	Limited internal volume constrains battery size and future addons
Manufacturability	Simpler tooling for assembly, but higher material cost	May require more precise tooling and inserts; lower material cost per unit
Weight	Higher overall weight (less mobile)	Lower weight improves portability and reduces shipping costs
Cable Management	Ample space for routing, strain relief, and connectors	Limited routing channels; requires careful cable planning and retention features

Table 1: Comparison between Initial and Final 3D enclosure designs.

3.5. IoT-Enabled Web Platform for IV Fluid Monitoring

To enhance the functionality, accessibility, and reliability of our IV fluid drop-rate monitoring device, we integrated it with an IoT-based web platform. This online system enables real-time monitoring of infusion status, patient information, and alert conditions from any location within the hospital network. The website acts as the central communication hub between the hardware device and the medical staff, ensuring continuous and accurate data visibility.

System Architecture Overview

The web platform is built using a modern, modular IoT architecture consisting of a Python backend, a React.js frontend, an SQL database, and REST-based API communication. Each component was selected to ensure reliability, scalability, and smooth real-time data handling.

3.5.1 Backend

The backend of the system is implemented using Python due to its strong ecosystem for IoT and data processing, high reliability for handling sensor data from the ESP module, seamless integration with databases and APIs, and fast development and debugging capabilities—crucial for rapid prototyping.

3.5.2 Frontend

The frontend website was developed using React.js, chosen for its ability to deliver fast, dynamic, real-time interfaces. It supports component-based development for modularity and maintainability. React's efficient state management ensures smooth live updates of drop rate, warnings, and patient data without requiring page reloads.



Purpose of the Web Platform




The purpose of developing a dedicated website for monitoring the IV fluid drop-rate device includes the following:



1. **Real-Time Monitoring:** Allows nurses to view the live drop rate and infusion status from any workstation without physically approaching the drip stand.
2. **Early Warning and Alerts:** If the drop rate deviates, slows, or stops, the website displays warning and critical alerts. This enables immediate medical action and reduces infusion errors.
3. **Remote Accessibility:** Staff can monitor multiple patients across different wards through one centralized dashboard.
4. **Data Logging and History:** All infusion sessions are stored, assisting with patient documentation, medical audits, and analysis of infusion trends and device performance.
5. **Improved Workflow Efficiency:** Reduces the need for manual drop counting, allowing medical staff to focus on more critical care activities.


6. Scalability: More devices can be added easily due to the API-based modular system design.



3.6. User Interface


Page Name	Image	Description
Load Page		The website is designed with a simple, clear, and secure interface to ensure smooth use by hospital staff. Each section of the interface supports a specific task in the IV fluid monitoring workflow, starting from logging in to monitoring patients and managing historical data.
Login Interface		The hospital administration begins by logging into the website through a secure login page. Password protection ensures that only authorized staff can access the system. This protects patient information, maintains privacy, and prevents unauthorized changes to infusion details or device settings. Security is essential in a hospital environment, and this interface ensures that all data is accessed only by verified users.

Page Name	Image	Description
Home Page Dashboard		The Home Page provides a quick overview of the system. It shows the total number of active assignments, total devices registered, total patients, and total beds in the system. From here, the admin can register a new device, register a patient, or register a bed to a specific ward. It also displays all ongoing active assignments, giving a clear snapshot of the current workload and device usage across the hospital.
Device Registration Interface		The Device Registration page is used when a new device is introduced to the system for the first time. The admin must enter the device's IP address and device number both of which are unique identifiers. This ensures that each device is mapped accurately and prevents conflicts when reading drop-rate data.
Patient Registration Interface		The Patient Registration interface allows the admin to enter and save patient details into the system. The required information includes the patient's full name, NIC number, and other relevant details. Registering the patient ensures that all infusion data is correctly linked to the correct individual in the hospital.

Page Name	Image	Description
Bed Registration Interface		<p>The Bed Registration page lets the admin assign a bed to a ward. A new bed can be added, or an existing bed can be reassigned to a different ward. The bed number and corresponding ward must be entered. This helps maintain an accurate map of where each patient and device is physically located.</p>
Administration Interface – Assignment Creation and Editing		<p>After all basic setups are complete, the Administration interface is used to create new infusion assignments. Here, a bed, a patient, and a device are linked together. The admin also enters the required drop rate, the drug name being infused, and the diagnosis. This page also allows editing an assignment which is useful when a previously assigned device needs to be moved to another patient or when patient details change.</p>

Page Name	Image	Description																														
All Devices Interface – Main Monitoring Page	 <table border="1"> <thead> <tr> <th>STATUS</th> <th>DEVICE</th> <th>IP ADDRESS</th> <th>PATIENT</th> <th>BED</th> <th>DRUG NAME</th> <th>EXPECTED RATE</th> <th>CURRENT RATE</th> <th>DEVIATION</th> <th>CONDITION</th> </tr> </thead> <tbody> <tr> <td>Offline</td> <td>ESP32_005</td> <td>172.28.10.4</td> <td>oshadha perera 200326111972</td> <td>BED_01 ICU</td> <td>alcohol</td> <td>17 drops/min</td> <td>.. drops/min</td> <td>.. drops/min</td> <td>Offline</td> </tr> <tr> <td>Offline</td> <td>ESP32_V3_001</td> <td>192.168.8.116</td> <td>Upani 200384100283</td> <td>BED_02 ICU</td> <td>N/A</td> <td>153 drops/min</td> <td>.. drops/min</td> <td>.. drops/min</td> <td>Offline</td> </tr> </tbody> </table> <p>Status Legend:</p> <ul style="list-style-type: none"> Normal: Within 5% of expected rate Warning: 5-15% deviation Critical: Above 15% deviation Offline: Device not responding 	STATUS	DEVICE	IP ADDRESS	PATIENT	BED	DRUG NAME	EXPECTED RATE	CURRENT RATE	DEVIATION	CONDITION	Offline	ESP32_005	172.28.10.4	oshadha perera 200326111972	BED_01 ICU	alcohol	17 drops/min	.. drops/min	.. drops/min	Offline	Offline	ESP32_V3_001	192.168.8.116	Upani 200384100283	BED_02 ICU	N/A	153 drops/min	.. drops/min	.. drops/min	Offline	<p>The All Devices page is the core of the system. It displays every device along with its current drop rate, required drop rate, and any deviations. The page shows whether a device is online or offline, its IP address, and the drug being administered. Warning and critical alerts are shown clearly, helping staff identify problems immediately. This is the most important page, as it provides the complete real-time monitoring overview for all active infusions.</p>
STATUS	DEVICE	IP ADDRESS	PATIENT	BED	DRUG NAME	EXPECTED RATE	CURRENT RATE	DEVIATION	CONDITION																							
Offline	ESP32_005	172.28.10.4	oshadha perera 200326111972	BED_01 ICU	alcohol	17 drops/min	.. drops/min	.. drops/min	Offline																							
Offline	ESP32_V3_001	192.168.8.116	Upani 200384100283	BED_02 ICU	N/A	153 drops/min	.. drops/min	.. drops/min	Offline																							

Page Name	Image	Description
Device Monitor – Individual Device View		<p>The Device Monitor interface allows users to search by device number, patient name, or bed number to view one specific ongoing assignment. It shows full patient details along with the live drip count and deviation indicators.</p> <p>Percentage = Deviation x 100% / required drop rate</p> <p>Normal (Green): Percentage \pm 5% Warning (Yellow): Percentage \pm 15% Critical (Red): Percentage \pm 15% for 5 consecutive readings</p> <p>At the bottom, a real-time graph displays the drop rate versus time. The required drop rate is shown as a straight green line, while the live drop rate is plotted dynamically, allowing staff to quickly understand infusion stability and trends.</p>
Settings (Manage) Interface		<p>The Settings page provides full administrative control. Staff can edit device information, edit patient details, or update bed records. It also includes the option to delete any unnecessary or outdated information, helping maintain a clean and accurate database.</p>

Page Name	Image	Description																																													
Patient History Interface	 <table border="1"> <thead> <tr> <th>DATE & TIME</th> <th>EXPECTED RATE</th> <th>ACTUAL RATE</th> <th>DEVIATION</th> <th>STATUS</th> </tr> </thead> <tbody> <tr><td>11/25/2025, 10:29:23 AM</td><td>557 drops/min</td><td>210.4078 drops/min</td><td>346.59 drops/min</td><td>Critical</td></tr> <tr><td>11/25/2025, 9:49:19 AM</td><td>557 drops/min</td><td>116.6938 drops/min</td><td>440.31 drops/min</td><td>Critical</td></tr> <tr><td>11/25/2025, 9:47:47 AM</td><td>12 drops/min</td><td>159.679 drops/min</td><td>147.68 drops/min</td><td>Critical</td></tr> <tr><td>11/25/2025, 9:47:43 AM</td><td>12 drops/min</td><td>159.679 drops/min</td><td>147.68 drops/min</td><td>Critical</td></tr> <tr><td>11/25/2025, 9:47:37 AM</td><td>12 drops/min</td><td>159.679 drops/min</td><td>147.68 drops/min</td><td>Critical</td></tr> <tr><td>11/25/2025, 9:47:33 AM</td><td>12 drops/min</td><td>159.679 drops/min</td><td>147.68 drops/min</td><td>Critical</td></tr> <tr><td>11/25/2025, 9:47:27 AM</td><td>12 drops/min</td><td>170.3347 drops/min</td><td>158.33 drops/min</td><td>Critical</td></tr> <tr><td>11/25/2025, 9:47:23 AM</td><td>12 drops/min</td><td>2.338062 drops/min</td><td>9.66 drops/min</td><td>Critical</td></tr> </tbody> </table>	DATE & TIME	EXPECTED RATE	ACTUAL RATE	DEVIATION	STATUS	11/25/2025, 10:29:23 AM	557 drops/min	210.4078 drops/min	346.59 drops/min	Critical	11/25/2025, 9:49:19 AM	557 drops/min	116.6938 drops/min	440.31 drops/min	Critical	11/25/2025, 9:47:47 AM	12 drops/min	159.679 drops/min	147.68 drops/min	Critical	11/25/2025, 9:47:43 AM	12 drops/min	159.679 drops/min	147.68 drops/min	Critical	11/25/2025, 9:47:37 AM	12 drops/min	159.679 drops/min	147.68 drops/min	Critical	11/25/2025, 9:47:33 AM	12 drops/min	159.679 drops/min	147.68 drops/min	Critical	11/25/2025, 9:47:27 AM	12 drops/min	170.3347 drops/min	158.33 drops/min	Critical	11/25/2025, 9:47:23 AM	12 drops/min	2.338062 drops/min	9.66 drops/min	Critical	<p>The Patient History interface allows staff to search for a specific patient and view all previous infusion records. The data is shown both as a graph and in a grid format, including recorded conditions such as normal, warning, or critical states. An important design choice is that the system stores data at 1 minute intervals. Since the device sends continuous real-time readings, saving every single reading would overload the database and cause storage issues. This optimized method keeps the system fast, efficient, and reliable while preserving meaningful historical trends.</p>
DATE & TIME	EXPECTED RATE	ACTUAL RATE	DEVIATION	STATUS																																											
11/25/2025, 10:29:23 AM	557 drops/min	210.4078 drops/min	346.59 drops/min	Critical																																											
11/25/2025, 9:49:19 AM	557 drops/min	116.6938 drops/min	440.31 drops/min	Critical																																											
11/25/2025, 9:47:47 AM	12 drops/min	159.679 drops/min	147.68 drops/min	Critical																																											
11/25/2025, 9:47:43 AM	12 drops/min	159.679 drops/min	147.68 drops/min	Critical																																											
11/25/2025, 9:47:37 AM	12 drops/min	159.679 drops/min	147.68 drops/min	Critical																																											
11/25/2025, 9:47:33 AM	12 drops/min	159.679 drops/min	147.68 drops/min	Critical																																											
11/25/2025, 9:47:27 AM	12 drops/min	170.3347 drops/min	158.33 drops/min	Critical																																											
11/25/2025, 9:47:23 AM	12 drops/min	2.338062 drops/min	9.66 drops/min	Critical																																											
Log Out Option		<p>Finally, a simple logout button allows users to securely exit the system. This helps maintain privacy and ensures that only authorized staff continue to have access.</p>																																													

3.7. Final Product and Product Dissection

A PCB was not printed and hardware components were soldered on to a dot board.




Figure 1: Product Dissection



Figure 2: Final Product

4. Accuracy Testing Results

Name of Test	Objective	Procedure	Pass / Fail Criteria	Status
Volumetric Chrono-metric Accuracy Test	Check if measured drop rate matches actual infusion rate	Set target volume and infusion time on the device. Measure drops with both device and manual reference (stopwatch & count). Repeat at low, medium, and high flow rates.	Pass: Device reading within $\pm 3\%$ or ± 2 drops/min of reference at all flow rates. Fail: Deviation beyond this limit in any test condition.	Pass
Photonic Sensitivity Validation	Verify IR sensor detects all drops correctly	Generate controlled drops using a calibrated IV set. Record counts with device and independent reference. Test under different lighting conditions.	Pass: more than 99% drops detected, less than 1% false positives. Fail: Any frequent misdetection, missed drops, or high false alarms.	Pass
Alarm Responsiveness Verification	Confirm alarms trigger for deviations	Set required drop rate on the device. Introduce slow flow (20%), fast flow (+20%), and blockage (0 drops). Record buzzer, LED, and web-app alarm response time.	Pass: Alarm activates immediately after the deviation. Fail: Alarm is delayed, missing, or activates during normal operation.	Pass
Long-Duration Stability Test	Ensure accuracy over extended operation	Run the device continuously for 24 hrs for a 500 ml saline volume infusion. Compare readings with reference every 6 hrs. Monitor for drift and error accumulation.	Pass: Cumulative drift less than $\pm 2\%$ over 24 hrs and device remains functional without resets. Fail: Drift greater than $\pm 2\%$, frequent errors, or device malfunctions.	Pass
Power Continuity and Recovery Test	Test battery backup and recovery after power loss	Switch from AC to battery power during normal operation. Measure backup run-time until battery depletion. Restore AC power and check for data retention.	Pass: Backup lasts for more than 8 hrs, data/settings remain intact, operation resumes within 10 s. Fail: Backup less than 8 hrs, data lost/corrupted, or system fails to resume properly.	In Progress
Data Transmission Integrity Test	Ensure correct web app data logging	Connect device to network and start continuous transmission. Introduce unstable conditions (packet loss, disconnects). Compare hospital database logs with device logs.	Pass: Data delay less than 10 s, recovery after disconnection, no missing records. Fail: Frequent packet loss, incorrect timestamps, or missing records.	Pending

User Interaction Validation	Confirm ease of use and prevention of input errors	Input volume/time and start infusion. Test functions: alarm mute, acknowledge, and calibration. Enter incorrect values to test error detection.	Pass: more than 90% of tasks completed correctly; input errors such as unrealistic volume are flagged and prevented. Fail: Frequent operator mistakes without system warning, unclear interface, or failure to accept correct input.	Pass
-----------------------------	--	---	---	------

5. Final Budget and Pricing Strategy

Final budget was calculated by adding the cost of PCB printing.

Expense	Unit price (Rs.)	Quantity	Amount (Rs.)
ESP32 WiFi Bluetooth Dev Board	1390	1	1390
3x4 Matrix Keypad Module	840	1	840
7 Segment 4 Bits Digital LED Display Module	345	1	345
OLED Display 0.96" I2C	690	1	690
IR Infrared Obstacle Avoidance Sensor Module	125	1	125
AC DC Power adopter 5 V	590	1	590
DC power socket plug jack connector	60	1	60
6x2mm Neodymium Magnet	30	4	120
Buzzer	50	1	50
TP4056 charging module	50	1	50
Switches	190	2	380
Miscellaneou (wires etc.)			150
Component Cost			4790
PCB print	2980	1	2980
Enclosure Print	3650	1	3650
TOTAL COST			11420

Figure 3: Budget for production of one item

Cost Per Unit (Bulk Production) +100 Units		
Expenses	Amount per Item	
Component Cost (30-50% decrement)- 40% off	4790 x 0.60	2874
Batch 3D printing (+100)	800-1400	1100
PCB Print(+100)	300 - 600	450
Manufacturing and Assembly cost (+100) (soldering, wiring etc.)	400- 600	500
Packaging and Logistics (Box, Cushioning, labelling etc.)	300 -500	400
Cost per unit		5324
Software Implementation - One time Software Licence Fee		1500
Marketing Sales & Operational overhead (10-20%) - 15%	5324 x 0.15	799
Profit Margin (15-40%) - 20% Profit	(5324 + 799) x 0.2	1225
Estimated Total Price per unit		8850
TOTAL UNIT PRICE = Rs. 8850 .00		

Figure 4: Budget for bulk production of 100 units

6. Possible Regulatory Pathway

Regulatory Pathway	Nature / Required Submissions
Device Classification	Likely classified as a FDA Class II device. The device is non-invasive and is required to perform as expected and will not cause injury or harm to the user.
Ethical Approval	<p>Objective: Needed to conduct initial clinical trials.</p> <p>Required Submissions: Submit documents to obtain ethical clearance from the University Ethics Review Committee, University of Moratuwa (https://uom.lk/uerc).</p>
Preclinical Testing	<p>Objective: Validate functionality, safety, and usability.</p> <p>Required Submissions:</p> <ul style="list-style-type: none"> • Performance Testing: Validate functionality • Electrical Safety and EMC: Test compliance with IEC 60601-1 and IEC 60601-1-2 • Usability: Conduct human factors studies per FDA guidelines
NMRA Approval	<p>Objective: To conduct clinical trials.</p> <p>Required Submissions: Submit documents to obtain clearance from the National Medicines Regulatory Authority (NMRA), the responsible body for the regulation and control of registration, licensing, manufacture, importation and all other aspects pertaining to medical devices in Sri Lanka to further continue clinical trials.</p>
Software and Cybersecurity	<p>Objective: Comply with IEC 62304 for software development and include testing results in the submission. Address cybersecurity risks with FDA's cybersecurity guidance, such as data encryption and risk mitigation.</p>
Premarket Submission (510(k))	<p>Required Submissions:</p> <ul style="list-style-type: none"> • Device description and intended use • Performance and safety test results • Risk analysis per ISO 14971 • Labeling
Manufacturing and Post-market Compliance	<p>Objective / Requirements:</p> <ul style="list-style-type: none"> • Implement a QMS compliant with 21 CFR Part 820 (Quality System Regulation) • Register the device and manufacturing facility with the FDA • Monitor post-market performance and report adverse events via the Medical Device Reporting (MDR) system

7. Conclusions and Future Improvements

7.0.1 Conclusion

The key uniqueness of our solution is adding an external display of the current infusion rate (drops per minute) to the simplest gravity-based infusion pumps commonly used in hospitals which is even a lacking feature in advanced pumps currently used in the hospital system as rate can change with patient movements. Our system provides alarms whenever the infusion rate deviates from the user-set target, ensuring patient safety without requiring expensive advanced pumps that can cost around Rs. 100 000.00. This hybrid approach delivers the core benefits of advanced infusion pumps—real-time monitoring and deviation alerts—in an affordable, user-friendly package suitable for everyday hospital use.

Advanced IV fluid regulation systems are not affordable to the local hospital system. At the moment, no one is providing a low-cost but reliable IV regulator made for the Sri Lankan context. This creates a clear space in the market that our solution can fill.

In conclusion, our team approached the development of this infusion pump with careful attention to both component selection and overall design. Every decision was made with the goal of meeting the exact functional requirements while ensuring reliability, safety, and usability. To achieve this, we consciously prioritized essential features, streamlining the device by omitting additional “smart” functions that, while interesting, were not critical to its core operation. This approach allowed us to create a solution that is not only cost-effective but also user-friendly and accessible, ensuring that it can serve its purpose efficiently in a real-world clinical setting.

7.0.2 Future Improvements

1. The system will be enhanced by enabling ESP32-to-ESP32 communication within each hospital ward. This will allow the creation of a local LED indicator network connecting all infusion setups in the ward. Each bed will have a dedicated LED indicator that reflects the current drip status — normal, deviating, or alarm — in real time.

This offline communication model will ensure that the monitoring and alert functions remain fully operational even if internet connectivity fails, making the system highly suitable for rural hospitals and low-resource environments where reliable internet access cannot be guaranteed.

2. We plan to keep the battery attachment and BMS as future improvements, while for now, we are focusing on implementing a power-saving mode to extend battery life. This mode will selectively turn off high-power components, such as the Wi-Fi module, when they are not needed, reducing energy consumption. By combining this with optimized microcontroller sleep states and minimal sensor activity during idle periods, we aim to maximize the efficiency of the current system while leaving room for future enhancements like advanced battery monitoring and management.
3. We are planning to add a software-based lock mechanism to the keypad to prevent patients from changing the drop rate, ensuring safety and consistent operation. Additionally, we are enhancing the system’s security in the web application by implementing two-step verification, providing an extra layer of protection for user

accounts and sensitive data.

For the web application,

4. One major improvement is reducing the amount of data the staff has to enter manually. In our current system, nurses and administrators must re-enter patient details, device information, drug data, and other fields specifically for IV infusion monitoring. But in a real hospital environment, this information is already recorded at the time of patient admission through the hospital's existing management software (HMS/EMR systems). Re-entering all this data again only for infusion monitoring is not practical and increases workload.
5. Another future enhancement is to make our web application work as an add-on module that can communicate directly with the hospital's existing software using APIs or standardized healthcare data exchange protocols. With this integration, once a patient is admitted in the hospital system, all relevant information can automatically sync to our device-assignment interface. This allows staff to assign devices quickly without typing the same details again. Alternatively, the infusion assignment feature itself can be embedded as a new function inside the hospital's existing software. This makes the system more realistic, efficient, and ready for real-world deployment in a busy clinical environment.

8. Task Allocation

Name	Index Number	Contribution
GUNATHUNGA U.A.	230218G	<ul style="list-style-type: none"> • Initial feasibility check of IR drop sensing method • Developing the IR sensing module with required initial functionality and its 3D Design • Content creation for project pitching and documentation
PERERA W.A.L.S.	230487D	<ul style="list-style-type: none"> • Initial feasibility check of image processing method for drop sensing • Component selection • Microcontroller programming • Soldering and assembling • Documentation
ABEYWARDHANA T.C.W.	230012U	<ul style="list-style-type: none"> • Initial feasibility check of capacitive method for drop sensing • Initial sketch of 3D design for part of the device fitting to the pole • 3D design development using Fusion • Documentation
PERERA K.W.A.O.V.	230481E	<ul style="list-style-type: none"> • Developing the web application • IoT integration of the device • Soldering and assembling • Documentation

9. Program for ESP32 Microcontroller

```

1 #include <Keypad.h>
2 #include <TM1637Display.h>
3 #include <Wire.h>
4 #include <Adafruit_SSD1306.h>
5 #include "esp_system.h"
6 #include <WiFi.h>
7 #include <WebServer.h>
8 #include <ArduinoJson.h>
9
10 const byte ROWS = 3;
11 const byte COLS = 4;
12
13 char keys[ROWS][COLS] = {
14     {'1','4','7','*'},
15     {'2','5','8','0'},
16     {'3','6','9','#'}
17 };
18
19 byte rowPins[ROWS] = {27, 26, 25};
20 byte colPins[COLS] = {33, 32, 14, 12};
21
22 Keypad keypad = Keypad(makeKeymap(keys), rowPins, colPins, ROWS,
23                         COLS);
24
25 #define CLK 2
26 #define DIO 15
27 TM1637Display display(CLK, DIO);
28
29 #define SCREEN_WIDTH 128
30 #define SCREEN_HEIGHT 64
31 #define OLED_RESET -1
32 #define SDA_PIN 21
33 #define SCL_PIN 22
34 #define BUZZER_PIN 13
35 Adafruit_SSD1306 oled(SCREEN_WIDTH, SCREEN_HEIGHT, &Wire,
36                         OLED_RESET);
37
38 // ===== WiFi CONFIGURATION =====
39 const char* ssid = "Dialog 4G - Indeepa"; // Change this to your
40 // WiFi name
41 const char* password = "Indeepa@05"; // Change this to your WiFi
42 // password
43 String deviceNumber = "ESP32_V3_001"; // Unique device identifier
44 String deviceIP = "";
45 WebServer server(80);
46
47 #define DROP_PIN 4
48 volatile unsigned long dropCount = 0;
49 volatile unsigned long lastDropTime = 0; // microseconds
50
51 bool patternPhase = false;
52 bool continuousPhase = false;
53
54 unsigned long melodyPrev = 0;
55 unsigned long panicPrev = 0;
56
57 bool inMelody = false;
58 bool inPanic = false;
59
60 int melodyStep = 0;
61 int panicStep = 0;
62
63 unsigned long medicineVolume = 0; // mL
64 unsigned long durationMinutes = 0;
65 float presetDropRate = 0; // DPM
66 unsigned long starttime = 0;
67 unsigned long alarmBlockUntil = 0; // Track start time for
       runtime calculations
68
69 // buzz blinking
70 unsigned long previousBlinkTime = 0;
71 bool buzzState = false;
72 const unsigned long blinkInterval = 250; // ms
73
74 // display updating
75 unsigned long previousDisplayMillis = 0;
76 const unsigned long displayInterval = 500;
77
78 // Beeping overtime correction
79 unsigned long overLimitStart = 0;
80 bool alarmActive = false;
81 const unsigned long overLimitDelay = 10000;
82
83 #define NUM_SAMPLES 5
84 volatile unsigned long dropTimes[NUM_SAMPLES];
85 volatile int dropIndex = 0;
86 bool enoughDrops = false;
87
88 float currentDPM = 0.0;
89
90 // ===== FUNCTION DECLARATIONS =====
91 void connectWiFi();
92 void setupWebServer();
93 void handleGetData();
94 void handleRoot();
95 void handleNotFound();
96
97 // ISR drop detection
98 void IRAM_ATTR onDropDetected() {
99     unsigned long now = micros();
100
101    if (now - lastDropTime > 150000) {
102        dropTimes[dropIndex] = now;
103        dropIndex = (dropIndex + 1) % NUM_SAMPLES;
104        dropCount++;
105        lastDropTime = now;
106
107        if (dropCount >= NUM_SAMPLES) enoughDrops = true; // start
               DPM calculation
108    }
109 }
110
111 // keyboard input
112 unsigned long getUserInput(String prompt) {
113     String input = "";
114     oled.clearDisplay();
115     oled.setTextSize(2);
116     oled.setTextColor(SSD1306_WHITE);
117     oled.setCursor(0,0);
118     oled.println(prompt);
119     oled.display();
120
121     while (true) {
122         server.handleClient(); // Handle web requests during input
123         char key = keypad.getKey();
124         if (key) {
125             if (key >= '0' && key <= '9') {
126                 input += key;
127             } else if (key == '*') { // backspace
128                 if (input.length() > 0) input.remove(input.length() -
1);
129             } else if (key == '#') { // confirm
130                 if (input.length() > 0) return input.toInt();
131             }
132
133             // Display current input
134             oled.clearDisplay();
135             oled.setTextSize(2);
136             oled.setCursor(0,0);
137             oled.println(prompt);
138             oled.setTextSize(3);
139             oled.setCursor(0,30);
140             oled.println(input);
141             oled.display();
142         }
143     }
144 }
145
146 // --- Calculate instantaneous DPM using last few drops ---
147 float calculateInstantDPM() {
148     if (!enoughDrops) return 0; // wait until 5 drops detected
149
150     // Find oldest and newest drop in buffer
151     unsigned long oldest = dropTimes[dropIndex]; // circular
               buffer
152     unsigned long newest = dropTimes[(dropIndex + NUM_SAMPLES - 1) %
NUM_SAMPLES];
153
154     unsigned long dt = newest - oldest; // total microseconds
               between 5 drops
155     if (dt == 0) return 0;
156
157     float dpm = (NUM_SAMPLES - 1) * 60.0 * 1000000.0 / dt; // extrapolate to drops per minute
158     return dpm;
159 }
160
161 void playMelodyPattern() {
162     static int sequence[] = {1,0,1,0,0,1,1,0}; // melody feel
163     static int len = 8;
164
165     if (millis() - melodyPrev > 200) { // slow beautiful
166         melodyPrev = millis();
167         melodyStep = (melodyStep + 1) % len;
168         digitalWrite(BUZZER_PIN, sequence[melodyStep]);
169     }
170 }
171
172 void playPanicPattern() {
173     const unsigned long interval = 80; // constant fast beeping
               without overload
174     if (millis() - panicPrev > interval) {
175         panicPrev = millis();
176         digitalWrite(BUZZER_PIN, !digitalRead(BUZZER_PIN));
177     }
178 }
179
180
181 void setup() {
182     Serial.begin(115200);
183     pinMode(BUZZER_PIN, OUTPUT);
184     digitalWrite(BUZZER_PIN, LOW);
185 }
```

```

186 // --- TM1637 ---
187 display.setBrightness(0x0f);
188 display.clear();
189
190 // --- OLED ---
191 Wire.begin(SDA_PIN, SCL_PIN);
192 if(!oled.begin(SSD1306_SWITCHCAPVCC, 0x3C)) {
193     Serial.println("OLED init failed");
194     for(;;);
195 }
196 oled.clearDisplay();
197 oled.display();
198 alarmBlockUntil = millis() + 10000; // 5 minutes
199
200 // Welcome screen
201 oled.clearDisplay();
202 oled.setTextSize(1);
203 oled.setTextColor(SSD1306_WHITE);
204 oled.setCursor(0,0);
205 oled.println("-----");
206 oled.setCursor(0,10);
207 oled.println(" Team InfuCare V3");
208 oled.setCursor(0,20);
209 oled.println(" IV Drop Detector");
210 oled.setCursor(0,30);
211 oled.println(" WiFi Version");
212 oled.setCursor(0,40);
213 oled.println("-----");
214 oled.display();
215 delay(2000);
216
217 // Connect to WiFi
218 oled.clearDisplay();
219 oled.setTextSize(1);
220 oled.setCursor(0,10);
221 oled.println("Connecting WiFi...");
222 oled.display();
223 connectWiFi();
224
225 if (WiFi.status() == WL_CONNECTED) {
226     deviceIP = WiFi.localIP().toString();
227     oled.setCursor(0,20);
228     oled.println("Connected!");
229     oled.setCursor(0,30);
230     oled.println(deviceIP.c_str());
231     oled.display();
232     delay(2000);
233 }
234 else {
235     oled.setCursor(0,20);
236     oled.println("WiFi Failed");
237     oled.setCursor(0,30);
238     oled.println("Offline Mode");
239     oled.display();
240     delay(2000);
241 }
242
243 // Setup web server
244 setupWebServer();
245
246 // --- Drop Sensor Interrupt ---
247 pinMode(DROP_PIN, INPUT_PULLUP); // active LOW
248 attachInterrupt(digitalPinToInterrupt(DROP_PIN),
                 onDropDetected, FALLING);
249
250 // --- Get medicine volume and duration from user ---
251 medicineVolume = getUserInput("Med Vol (mL):");
252 durationMinutes = getUserInput("Duration (min):");
253
254 // Calculate preset drop rate: drops per minute
255 presetDropRate = (medicineVolume * 20.0) / durationMinutes;
256
257 oled.clearDisplay();
258 oled.setTextSize(2);
259 oled.setCursor(0,0);
260 oled.println("Preset DPM");
261 oled.setTextSize(3);
262 oled.setCursor(0,30);
263 oled.println(presetDropRate);
264 oled.display();
265
266 delay(3000); // show preset rate for 3 seconds
267
268 starttime = millis(); // Start timing
269
270 Serial.println("\n=====");
271 Serial.println("Device Ready!");
272 Serial.println("Device Number: " + deviceNumber);
273 Serial.println("IP Address: " + deviceIP);
274 Serial.println("Expected Drop Rate: " + String(presetDropRate,
275                                         1) + " drops/min");
276 Serial.println("Data Endpoint: http://" + deviceIP + "/data");
277 }
278
279 void loop() {
280     // --- Handle web server requests ---
281     server.handleClient();
282 }

283 // --- Keypad handling ---
284 char key = keypad.getKey();
285 if (key == '*') {
286     oled.clearDisplay();
287     oled.setTextSize(2);
288     oled.setCursor(0,0);
289     oled.println("Rebooting...");
290     oled.display();
291     digitalWrite(BUZZER_PIN, LOW);
292     delay(500);
293     esp_restart();
294 }
295 if (key == '#') {
296     dropCount = 0; // reset counter
297     display.clear();
298     oled.clearDisplay();
299     oled.setTextSize(2);
300     oled.setCursor(0,0);
301     oled.println("Reset Done");
302     digitalWrite(BUZZER_PIN, LOW);
303     oled.display();
304     delay(100);
305 }
306
307 unsigned long currentMillis = millis();
308
309 // --- Calculate instantaneous current DPM ---
310 currentDPM = calculateInstantDPM();
311 float percentageError = 0;
312
313 if(presetDropRate > 0){
314     percentageError = fabs(currentDPM - presetDropRate) /
315                         presetDropRate * 100.0;
316 }
317
318 // Track when deviation starts
319 static bool deviationStarted = false;
320
321 if (percentageError >= 10) {
322     if (!deviationStarted) {
323         overLimitStart = millis(); // mark start of deviation
324         deviationStarted = true;
325     }
326 } else {
327     deviationStarted = false; // back to normal
328 }
329
330 if(millis() < alarmBlockUntil){
331     digitalWrite(BUZZER_PIN, LOW);
332     oled.clearDisplay();
333     oled.setTextSize(2);
334     oled.setCursor(0,44);
335     oled.println("ADJUSTING");
336     oled.display();
337 }
338
339 // --- Update displays every 0.5 second ---
340 if (currentMillis - previousDisplayMillis >= displayInterval)
341 {
342     // TM1637: current DPM
343     display.showNumberDec((int)currentDPM, false);
344
345     // OLED: preset DPM, current DPM, and WiFi status
346     oled.clearDisplay();
347     oled.setTextSize(1);
348     oled.setCursor(0,0);
349     oled.println("Preset DPM:");
350     oled.setTextSize(3);
351     oled.setCursor(0,18);
352     oled.println((int)presetDropRate);
353     oled.setTextSize(1);
354     oled.setCursor(0,40);
355     oled.println("Current DPM:");
356     oled.setTextSize(2);
357     oled.setCursor(0,50);
358     oled.println((int)currentDPM);
359
360     // WiFi status at bottom
361     oled.setTextSize(1);
362     oled.setCursor(0,54);
363     if(WiFi.status() == WL_CONNECTED) {
364         // oled.println("WiFi: " + deviceIP.substring(deviceIP.
365                     lastIndexOf('.') + 1));
366     } else {
367         oled.println("WiFi: Off");
368     }
369     oled.display();
370
371     previousDisplayMillis = currentMillis;
372 }
373
374 // not enough drops do nothing
375 if(!enoughDrops){
376     digitalWrite(BUZZER_PIN, LOW);
377     return;
378 }

```

```

379     if (percentageError < 10) {
380         // NORMAL
381         digitalWrite(BUZZER_PIN, LOW);
382     } else if (percentageError < 20) {
383         // WARNING
384         if (millis() - melodyPrev > 500) {
385             melodyPrev = millis();
386             digitalWrite(BUZZER_PIN, !digitalRead(BUZZER_PIN));
387         }
388     } else {
389         // CRITICAL
390         unsigned long devTime = millis() - overLimitStart;
391
392         if (devTime >= deviationDelay + patternDuration) {
393             playPanicPattern(); // panic
394         } else if (devTime >= deviationDelay) {
395             playMelodyPattern(); // melody stage
396         } else {
397             digitalWrite(BUZZER_PIN, LOW); // initial deviation
398         }
399     }
400 }
401
402 // ===== WiFi CONNECTION =====
403 void connectWiFi() {
404     Serial.print("Connecting to WiFi: ");
405     Serial.println(ssid);
406
407     // Ensure station mode before connecting
408     WiFi.mode(WIFI_STA);
409     WiFi.begin(ssid, password);
410
411     int attempts = 0;
412     while (WiFi.status() != WL_CONNECTED && attempts < 30) {
413         delay(500);
414         Serial.print(".");
415         attempts++;
416     }
417
418     if (WiFi.status() == WL_CONNECTED) {
419         Serial.println("\n WiFi Connected!");
420         Serial.print("IP Address: ");
421         Serial.println(WiFi.localIP());
422     } else {
423         Serial.println("\n WiFi Connection Failed!");
424         Serial.println("Device will work in offline mode.");
425     }
426 }
427
428 // ===== WEB SERVER SETUP =====
429 void setupWebServer() {
430     server.on("/", handleRoot);
431     server.on("/data", HTTP_GET, handleGetData);
432     server.onNotFound(handleNotFound);
433     server.begin();
434     Serial.println(" Web server started on port 80");
435 }
436
437 // ===== WEB SERVER HANDLERS =====
438 void handleRoot() {
439     String html = "<!DOCTYPE html><html><head><title>ESP32 V3 IV
        Monitor</title>";
440     html += "<meta name='viewport' content='width=device-width,
        initial-scale=1'>";
441     html += "<style>body{font-family:Arial;padding:20px;background:#f0f0f0;margin:0;}";
442     html += ".container{background:white;padding:20px;border-radius:
        8px;max-width:600px;margin:0 auto;box-shadow:0 2px 4px
        rgba(0,0,0,0.1)}";
443     html += "h1{color:#203f76;margin-top:0;}";
444     html += "table{width:100%;border-collapse:collapse;margin:20px
        0;}";
445     html += "td{padding:12px 8px;border-bottom:1px solid #ddd;}";
446     html += "td:first-child{font-weight:bold;color:#203f76;width
        :40%;}";
447     html += ".status{display:inline-block;padding:4px 12px;border-
        radius:12px;font-size:14px;}";
448     html += ".online{background:#d4edda;color:#155724;}";
449     html += ".offline{background:#ff8d7d;color:#721c24;}";
450     html += ".warning{background:#fff3cd;color:#856404;}";
451     html += "a{color:#57bdd4;text-decoration:none;font-weight:bold;}";
452     html += "a:hover{text-decoration:underline;}</style></head>";
453     html += "<body><div class='container'>";
454     html += "<h1> InfuCare IV Monitor V3</h1>";
455     html += "<table>";
456     html += "<tr><td>Device Number:</td><td>" + deviceNumber + "</td
        ></tr>";
457     html += "<tr><td>IP Address:</td><td>" + deviceIP + "</td></tr>";
458     html += "<tr><td>Status:</td><td><span class='status online'>
        Active</span></td></tr>";
459     html += "<tr><td>WiFi:</td><td><span class='status '" + String(
        WiFi.status() == WL_CONNECTED ? "online" : "offline") + " Disconnected"
        </span></td></tr>";
460
461     html += "<tr><td>Expected Rate:</td><td><strong>" + String(
        presetDropRate, 1) + "</strong> drops/min</td></tr>";
462     html += "<tr><td>Current Rate:</td><td><strong>" + String(
        currentDPM, 1) + "</strong> drops/min</td></tr>";
463     html += "<tr><td>Total Drops:</td><td><strong>" + String(
        dropCount) + "</strong></td></tr>";
464     html += "<tr><td>Volume Setting:</td><td>" + String(
        medicineVolume) + " ml</td></tr>";
465     html += "<tr><td>Time Setting:</td><td>" + String(
        durationMinutes) + " minutes</td></tr>";
466
467     float elapsedSeconds = (millis() - starttime) / 1000.0;
468     int minutes = int(elapsedSeconds / 60);
469     int seconds = int(elapsedSeconds) % 60;
470     html += "<tr><td>Running Time:</td><td>" + String(minutes) + "m
        " + String(seconds) + "s</td></tr>";
471
472     float deviation = abs(currentDPM - presetDropRate);
473     float deviationPercent = (presetDropRate > 0) ? (deviation /
        presetDropRate) * 100 : 0;
474     html += "<tr><td>Deviation:</td><td>" + String(deviation, 1) + "
        drops/min (" + String(deviationPercent, 1) + "%)</td></tr>";
475
476     // Alarm status
477     String alarmStatus = "Normal";
478     if(currentDPM > presetDropRate) {
479         unsigned long devTime = millis() - overLimitStart;
480         if(devTime >= deviationDelay + patternDuration) {
481             alarmStatus = "PANIC ALARM";
482         } else if(devTime >= deviationDelay) {
483             alarmStatus = "Melody Alert";
484         } else {
485             alarmStatus = "Deviation Detected";
486         }
487     }
488     html += "<tr><td>Alarm Status:</td><td><span class='status '" +
        String(alarmStatus == "Normal" ? "online" : "warning") + "</span></td></tr>";
489
490     html += "</table>";
491     html += "<p style='margin-top:20px;text-align:center;'><a href
        ='#data'> View JSON Data</a></p>";
492     html += "<p style='text-align:center;color:#666;font-size:12px;
        margin-top:30px;'>InfuCare Team 2025 - V3</p>";
493     html += "</div></body></html>";
494
495     server.send(200, "text/html", html);
496     Serial.println(" Web page served");
497 }
498
499 void handleGetData() {
500     // Prepare JSON response
501     StaticJsonDocument<512> doc;
502     doc["device_number"] = deviceNumber;
503     doc["drop_rate"] = currentDPM;
504     doc["drop_count"] = dropCount;
505     doc["expected_drop_rate"] = presetDropRate;
506     doc["timestamp"] = String(millis());
507     doc["device_status"] = "online";
508     doc["wifi_status"] = (WiFi.status() == WL_CONNECTED) ? "connected" :
        "disconnected";
509     doc["ip_address"] = deviceIP;
510
511     // Calculate deviation
512     float deviation = abs(currentDPM - presetDropRate);
513     float deviationPercent = (presetDropRate > 0) ? (deviation /
        presetDropRate) * 100 : 0;
514     doc["deviation"] = deviation;
515     doc["deviation_percent"] = deviationPercent;
516
517     // Add settings
518     doc["volume_setting"] = medicineVolume;
519     doc["time_setting"] = durationMinutes;
520
521     // Add runtime
522     float elapsedSeconds = (millis() - starttime) / 1000.0;
523     doc["runtime_seconds"] = elapsedSeconds;
524
525     // Alarm status
526     String alarmStatus = "normal";
527     if(currentDPM > presetDropRate) {
528         unsigned long devTime = millis() - overLimitStart;
529         if(devTime >= deviationDelay + patternDuration) {
530             alarmStatus = "panic";
531         } else if(devTime >= deviationDelay) {
532             alarmStatus = "melody";
533         } else {
534             alarmStatus = "deviation";
535         }
536     }
537     doc["alarm_status"] = alarmStatus;
538
539     String jsonData;
540     serializeJson(doc, jsonData);
541
542     // Send CORS headers
543     server.sendHeader("Access-Control-Allow-Origin", "*");
544
545 }
```

```
548     server.sendHeader("Access-Control-Allow-Methods", "GET");
549     server.send(200, "application/json", jsonData);
550
551     Serial.println(" Data request served - Rate: " + String(
552         currentDPM, 1) + " d/m");
553
554     void handleNotFound() {
555         server.send(404, "text/plain", "404: Not Found");
556     }
```

References

- [1] P. G. Yock, S. Zenios, J. Makower, T. J. Brinton, U. N. Kumar, F. T. J. Watkins, L. Denend, T. M. Krummel, and C. Q. Kurihara, Eds., Biodesign: The Process of Innovating Medical Technologies, 2nd ed. Cambridge, U.K.: Cambridge Univ. Press, 2015
- [2] O. Hoffman and O. Bacon, “Infusion Pumps,” in Making Healthcare Safer III: A Critical Analysis of Existing and Emerging Patient Safety Practices, K. K. Hall, S. Shoemaker-Hunt, L. Hoffman, et al., Eds. Rockville, MD, USA: Agency for Healthcare Research and Quality (US), Mar. 2020. [Online]. Available: <https://www.ncbi.nlm.nih.gov/books/NBK555506/>
- [3] B. Jung, K. S. Seo, S. J. Kwon, K. Lee, S. Hong, H. Seo, G. Y. Kim, G. M. Park, J. Jeong, and S. Seo, “Efficacy evaluation of syringe pump developed for continuous drug infusion,” J. Dent. Anesth. Pain Med., vol. 16, no. 4, pp. 303–307, Dec. 2016. doi: 10.17245/jdapm.2016.16.4.303. PMID: 28879319; PMCID: PMC5564196. [Online]. Available: <https://pmc.ncbi.nlm.nih.gov/articles/PMC5564196/>